



# Flexidose light

## The dosimetry management solution

- **FlexiDose Light** is a dosimetry management software allowing the assignment of the dosimeters, the dose plurality and the workers history. **Flexidose Light** is perfectly adapted for small sets of dosimeters. Designed for small organizations on a customer-server basis, it can control up to 10 dosimeter readers.



### ➤ Technical characteristics

#### Management / Performances

- Up to 10 readers
- 2000 movements per day
- Up to 100 people
- Compatible with Saphydose γ1 dosimeter
- Compatible with LMF readers (option) and SAPHYR reader

#### Functionalities

- Data access managed by identifier/password (Two profiles: system administrator and user)
- Real time access control, an access authorization can be regulated individually, by id, work or period
- Report and export integrated tools transfer data to another system (reports available in XML or CSV)

#### Environment :

- SQL data base (Access compatible). Operates on Windows 2000 or more recent
- The same PC can integrate the communication with the dosimeter, the data base and the software.

- **User-friendly Interface**
- **Data access managed by identifier / password**
- **Real-time access control**
- **Automatic detection of the readers or of network problems**
- **Integrated reports and export tools**
- **SQL Server standard database**
- **Management of externalized agents.**